Tube Guide



A full line of BD Vacutainer Blood Collection Needles, Needle Holders and Blood Collection Sets is also available.

BD Vacutaine Tubes With Hemogard™ C		BD Vacutainer™ Tubes With Conventional Stopper		nversions at Blood Collection*	Laboratory Use	Your Lab's Draw Volume/Remark
	Gold	Red/ Black	Clot activator and gel for serum separation	5	BD Vacutainer™ SST™ Tube for serum determinations in chemistry. Tube inversions ensure mixing of clot activator with blood. Blood clotting time 30 minutes.	
	Light Green	Green/ Gray	Lithium heparin and gel for plasma separation	8	BD Vacutainer™ PST™ Tube for plasma determinations in chemistry. Tube inversions prevent clotting.	
	Red		 None (glass) Clot activator (plastic tube with Hemogard closure) 	0 5	For serum determinations in chemistry and serology. Glass serum tubes are recommended blood banking. Plastic tubes contain clot activate and are not recommended for blood banking. Tube inversions ensure mixing of clot activator with blood and clotting within 60 minutes.	for ror
	Orange	Gray/ Yellow	• Thrombin	8	For stat serum determinations in chemistry. Tubinversions ensure complete clotting which usua occurs in less than 5 minutes.	oe Ally
	Royal Blue		 Sodium heparin Na₂EDTA None (serum tube) 	8 8 0	For trace-element, toxicology and nutritional-chemistry determinations. Special stopper formulation provides low levels of trace elements (see package insert).	
	Green		Sodium heparinLithium heparin	8 8	For plasma determinations in chemistry. Tube inversions prevent clotting.	
	Gray		 Potassium oxalate/ sodium fluoride Sodium fluoride/Na₂ EDTA Sodium fluoride (serum tub 	8 8 e) 8	For glucose determinations. Oxalate and EDTA anticoagulants will give plasma samples. Sodium fluoride is the antiglycolytic agent. Tube inversions ensure proper mixing of additional blood.	<i>r</i> e
	Tan		 Sodium heparin (glass) K₂EDTA (plastic) 	8 8	For lead determinations. This tube is certified to contain less than .01 µg/mL(ppm) lead. Tube inversions prevent clotting.	0
			Sodium polyapothal sulfapata (SDS)	8	SPS for blood culture specimen collections in	
	Yellow		polyanethol sulfonate (SPS) • Acid citrate dextrose additives (ACD): Solution A - 22.0g/L trisodium citrate, 8.0g/L citric acid, 24.5g/L dextrose Solution B - 13.2g/L trisodium citrate, 4.8g/L citric acid, 14.7g/L dextrose	8	micro-biology. Tube inversions prevent clotting ACD for use in blood bank studies, HLA phenotyping, DNA and paternity testing.	
	Lavender		 Liquid K₃EDTA (glass) Spray-dried K₂EDTA (plastic) 	8 8	$\rm K_3EDTA$ for whole blood hematology determin tions. $\rm K_2EDTA$ for whole blood hematology determinations and immunohematology testin (ABO grouping, Rh typing, antibody screening) Tube inversions prevent clotting.	g
	Pink		• Spray-dried K ₂ EDTA	8	For whole blood hematology determinations a immunohematology testing (ABO grouping, Rt typing, antibody screening). Designed with specross-match label for required patient informat by the AABB. Tube inversions prevent clotting.	n ecial
	Light Blue		 105M sodium citrate (≈3.2% 129M sodium citrate (3.8%) Citrate, theophylline, adenosine, dipyridamole (CTAD) 	3-4 3-4 3-4	For coagulation determinations. NOTE: Certain tests may require chilled specimens. Follow you institution's recommended procedures for collection and transport. CTAD for selected platelet function assays and routine coagulatior determination. Tube inversions prevent clotting	r 1
rtial-draw nl and 3 ml: 1		Small-volume Pediatric Tul (2ml: 10.25 x 47 mm, 3ml: 10.25 x	pes			
	Red		• None	0	For serum determinations in chemistry and serce Glass serum tubes are recommended for blood banking. Plastic tubes contain clot activator and not recommended for blood banking. Tube inversions ensure mixing of clot activator with bland clotting within 60 minutes.	I are
	Green	8	Sodium heparinLithium heparin	8 8	For plasma determinations in chemistry. Tube inversions prevent clotting.	
	Lavender	8	 Liquid K₃EDTA (glass) Spray-dried K₂EDTA (plastic) 	8 8	K ₃ EDTA for whole blood hematology determin tions. K ₂ EDTA for whole blood hematology determinations and immunohematology testin (ABO grouping, Rh typing, antibody screening) Tube inversions prevent clotting.	q
	Light Blue		 .105M sodium citrate (≈3.2%) .129M sodium citrate (3.8%) 	3-4	For coagulation determinations. Tube inversion prevent clotting. NOTE: Certain tests may required chilled specimens. Follow your institution's recommended procedures for collection and	s ire

www.bd.com